
























IMPLANT TORQUE SPECIFICATION GUIDE

MANUFACTURER	IMPLANT	TORQUE (NCM)	TORQUE SETTINGS	DRIVER & HEAD
ASTRA				
	ASTRA 3.0	15 NCM	1.27MM HEX	
	ASTRA 3.5-4.0	20 NCM	1.27 MM HEX	
	ASTRA 4.5-5.0	25 NCM	1.27MM HEX	
	ASTRA EV 3.0, 3.6, 4.2, 4.8, 5.4	25 NCM	1.27 MM HEX	
BIOHORIZONS				
	INTERNAL 3.5, 4.5, 5.7	30 NCM	1.27MM HEX	
	TAPERED INTERNAL 3.5, 4.5, 5.7	30 NCM	1.27MM HEX	
	EXTERNAL IMPLANT 3.5	30 NCM	1.27MM HEX	
	EXTERNAL IMPLANT 4.0, 5.0	30 NCM	1.27MM HEX	
NOBEL				
	ASC	35 NCM	NOBEL OMNIGRIP	
	ACTIVE 3.0	15 NCM	NOBEL UNIGRIP	
	ACTIVE/REPLACE CONICAL CONNECTION NP 3.5	35 NCM	NOBEL UNIGRIP	
	ACTIVE/REPLACE CONICAL CONNECTION RP 4.3, 5.0	35 NCM	NOBEL UNIGRIP	
	REPLACE NP 3.5	35 NCM	NOBEL UNIGRIP	
	REPLACE RP 4.3	35 NCM	NOBEL UNIGRIP	
	REPLACE WP 5.0	35 NCM	NOBEL UNIGRIP	
	REPLACE 6.0	35 NCM	NOBEL UNIGRIP	
	BRANEMARK NP 3.3	35 NCM	NOBEL UNIGRIP	
	BRANEMARK RP 3.75, 4.0	35 NCM	NOBEL UNIGRIP	
	BRANEMARK WP 5.0, 5.5	35 NCM	NOBEL UNIGRIP	
	MUA	15 NCM	NOBEL UNIGRIP	
STRAUMANN				
	BONE LEVEL 3.3 NC	35 NCM	SCS SCREWDRIVER	



IMPLANT TORQUE SPECIFICATION GUIDE

MANUFACTURER	IMPLANT	TORQUE (NCM)	TORQUE SETTINGS	DRIVER & HEAD
<i>STRAUMANN</i>				
	BONE LEVEL 4.1, 4.8 RC	35 NCM	SCS SCREWDRIVER	★
	STANDARD/STANDARD PLUS 4.8 RN	35 NCM	SCS SCREWDRIVER	★
	STANDARD/STANDARD PLUS 4.8 WN	35 NCM	SCS SCREWDRIVER	★
	STANDARD PLUS 3.3 NN	35 NCM	SCS SCREWDRIVER	★
	TISSUE LEVEL	30 NCM	STAR/TORX	
	MUA	15 NCM	SCS SCREWDRIVER	★
<i>DENTSPLY</i>				
	OSSEOSPEED EV 3.0 GREEN	25 NCM	0.050" HEX	⬡
	OSSEOSPEED EV 3.6 PURPLE	25 NCM	0.050" HEX	⬡
	OSSEOSPEED EV 4.2 YELLOW	25 NCM	0.050" HEX	⬡
	OSSEOSPEED EV 4.8 BLUE	25 NCM	0.050" HEX	⬡
	OSSEOSPEED EV 5.4 BROWN	25 NCM	0.050" HEX	⬡
	OSSEOSPEED PROFILE EV 4.2 YELLOW	25 NCM	0.050" HEX	⬡
	OSSEOSPEED PROFILE EV 4.8 BLUE	25 NCM	0.050" HEX	⬡
<i>GLIDEWELL</i>				
	INCLUSIVE TAPERED 3.7, 4.7, 5.2	35 NCM	0.050" HEX	⬡